September 5, 2013

Edan Medical

I had a two week period to evaluate your U50 ultrasound machine and, for the price point, I thought the over-all operation was quite good. I have listed my comments about this compact ultrasound machine below.

- 1) The Color Doppler was very good in sensitivity and spatial resolution for both superficial and deep vascular structures.
- 2) Very user- friendly system controls and functionality mode controls that are intuitive and easy to use.
- 3) B-Mode resolution is good for the curved-array transducer for abdominal and obstetrical scanning
- 4) This unit is small and light weight for portable and bedside examinations.

I feel that this system would be very competitive in any point-of care ultrasound lab.

Michael S. Cronan, RDMS

Chief Sonographer

U. C. Davis Medical Center

Sacramento, California